

1. **Part Number:** 06286-45
2. **Date Issued:** 09/07/22



3. Product Information:

Description:	Powder Lance, HDPE, 600mm length, 125mL; 20/PK
Use:	Single use, disposable powder sampler
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	601mm
Blade Length:	500mm
Diameter of Blade:	22mm
Nominal Weight of Sampler:	80g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

4. **Country of Origin:** Made and packed in the UK.
5. **Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.