

**1. Part Number:** 06286-47

**2. Date Issued:** 09/07/22



Description	Long Handled Spoon, HDPE, 10mL; 100/PK
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> </ul>
Method of Construction:	Single piece, injection moulded
Bowl Volume:	10ml
Handle Length:	285mm
Bowl Length:	77mm
Bowl Width:	29mm
Nominal Weight:	20g
Moulding & Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	100
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature