

1. Part Number: 06286-50

2. Photograph:



3. Product Information:

Description:	Scraper, Sterile, HDPE, 75mm width; 100/PK
Use:	Single use, disposable scraper
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Overall Length:	190mm
Blade Length:	75mm
Max. Blade Width:	75mm
Nominal Weight:	36g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (Minimum 25 kGy)
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

- 4. Country of Origin:** Made and packed in the UK.
- 5. Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.