

1. Part Number: 06286-52

2. Date Issued: 09/07/22



3. Product Information:

Description: Scraper, Sterile, HDPE, 110mm width; 100/PK

Use: Single use, disposable scraper

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Single piece, injection moulded

Overall Length: 225mm

Blade Length: 110mm

Blade Width: 110mm

Nominal Weight: 50g

Moulding & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (Minimum 25 kGy)

Quantity per Box: 100

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This device has not been ATEX tested

4. Country of Origin: Made and packed in the UK.

5. Repeated Use: This is a single-use, disposable sampler. Repeated use will lead to failure.