

# Nordson Divergent™ 1-Way Stopcocks Polycarbonate

Cole-Parmer®



Add a stopcock to your system for easy on/off control of your fluid. One-way stopcocks feature one female Luer fitting and a male Luer fitting with or without a lock. Luer fittings on stopcocks ensure a quick, leak-free connection. Diverter arrows indicate flow direction, enhancing operational efficiency and minimizing errors.



## CERTIFICATION

ISO 80369-7 compliant

**WARNING:** This product is not approved or intended for, and should not be used for, medical, clinical, surgical, or other patient-oriented applications.

## APPLICATIONS

- Chromatography, sampling ports, and single-stream titration
- Ideal for general-purpose, laboratory, and single-use applications

## FEATURES/BENEFITS

- Low-torque handle provides smooth adjustments
- Diverter arrows indicate flow direction
- DEHP free (non-DEHP)
- Transparent polycarbonate body makes it easy to see the flow

USA	+1.800.323.4340
	+1.847.549.7600
Canada	+1.800.363.5900
China	86.21.5109.9909
India	+1.800.266.1244

UK	+44 (0) 1480.272279
	+33 (0) 1.87.170142*
	+49 (0) 937.792030*
Italy	+39.02.84349215
All others	+1.847.549.7600

Cole-Parmer®

[coleparmer.com](http://coleparmer.com)

\*Inquiries from Germany and France are now handled in our St. Neots office by native-speaking experts.

### Specification and Ordering Table

Number of Ports	Port 1	Port 2	Max Operating Pressure (psi)	Item Number	Pack Size
2	Female Luer	Male Slip Luer	50	<a href="#">56140-96</a>	25
2	Female Luer	Male Lock Luer	50	<a href="#">56140-97</a>	25
2	Female Luer	Male Lock Luer	500	<a href="#">56140-99</a>	25

### FLOW PATH CONFIGURATIONS

