

# Cole-Parmer® 1-Way Stopcocks

## Natural Kynar® PVDF

Cole-Parmer®



Add a stopcock to your system for easy on/off control of your fluid. The handle provides a quick visual check of whether the valve is open or closed. Stopcocks feature luer fittings, making it easy to create a quick leak-free connection. Natural Kynar PVDF stopcocks are animal-derivative-, DEHP-, and BPA-free, making them ideal for bioprocessing research, pharmaceutical, and food and beverage applications.

### APPLICATIONS

- Perfect for titrations, chromatography, and cooling circuits
- Ideal for general-purpose, laboratory, and single-use applications

### FEATURES/BENEFITS

- Stopcocks provide accurate, leak-free flow control
- Handles rotate a full 360°
- Sterilize via autoclave, EtO, or gamma irradiation
- Exhibits exceptional chemical resistance and high purity
- Conforms to ISO 594-1 and 594-2 standards
- Free technical application support available to assist with product selection



### CERTIFICATION

USP Class VI, NSF 51, NSF 61, FDA, and RoHS compliant

USA	+1.800.323.4340
	+1.847.549.7600
Canada	+1.800.363.5900
China	86.21.5109.9909
India	+1.800.266.1244

UK	+44 (0) 1480.272279
	+33 (0) 1.87.170142*
	+49 (0) 937.792030*
Italy	+39.02.84349215
All others	+1.847.549.7600

Cole-Parmer®

[coleparmer.com](http://coleparmer.com)

\*Inquiries from Germany and France are now handled in our St. Neots office by native-speaking experts.

## Specification and Ordering Table

Port 1	Port 2	Material	Item Number	Pack Size
Female Luer	Female Luer	PVDF	<a href="#">50140-54</a>	10
Male Lock Luer	Female Luer	PVDF	<a href="#">50140-55</a>	10
Male Lock Luer	Male Lock Luer	PVDF	<a href="#">50140-56</a>	10
Male Lock Luer	Male Slip Luer	PVDF	<a href="#">50140-57</a>	10
Male Slip Luer	Female Luer	PVDF	<a href="#">50140-58</a>	10
Male Slip Luer	Male Slip Luer	PVDF	<a href="#">50140-59</a>	10