



## MB-800 Series Moisture Determination Balances

- Three-year industry leading warranty
- Readability as precise as 0.1 mg
- Touch screen display
- Automatic weighing chamber door



# Cole-Parmer® MB-800 Series Moisture Determination Balances

The MB-800 series of moisture balances combines state of the art features, an intuitive touch screen display, premium build quality, and Cole-Parmer's industry leading three-year warranty into an economical package that is an essential addition to any laboratory.

Onboard database storage and processing, and unparalleled connectivity in the form of RS-232, USB-A, USB-B, and Wi-Fi make the MB-800 series moisture determination balances the perfect choice for all your moisture analysis needs.

## Technical Specifications

Specification	MB-800-50.H	MB-800-50	MB-800-210
Max Capacity	50 g	50 g	210 g
Readability	0.1 mg	1 mg	1 mg
Calibration Type	External	External	External
Weigh Pan Dimensions	ø90 mm	ø90 mm	ø90 mm
Heater Type	Infrared	Infrared	Infrared
Max Temperature	160 °C	160 °C	160 °C
Device Dimensions	13 x 8 x 7.8 in (338 x 206 x 197 mm)		
Connectivity	2x RS-232, 2x USB-A (interchangeable), USB-B, Wi-Fi		

## Ordering Information

Description	Series No.	Model No.	Ordering No.
50 g, 0.1 mg Moisture Determination Balance with Touch Screen, High Resolution	MB-800	MB-800-50.H	10100-42
50 g, 1 mg Moisture Determination Balance with Touch Screen	MB-800	MB-800-50	10100-43
210 g, 1 mg Moisture Determination Balance with Touch Screen	MB-800	MB-800-210	10100-44

USA:  
+1.800.323.4340  
+1.847.549.7600

Canada: +1.800.363.5900  
China: +86.21.5109.9909  
France: +33 (0) 1486 37800  
Germany: +49 (0) 9377 92030

India: +1.800.266.1244  
Italy: +39 (0) 1313 89513  
UK: +44 (0) 1480 272279  
All other countries: +1.847.549.7600

**Cole-Parmer**