



Scales feature weighing, parts counting, percentage, checkweighing, dynamic, and accumulation functions.

***Withstands splashes,  
spills, and sprays***



## Features

- IP66-rated housing and indicator are suitable for washdown environments
- Rugged stainless-steel construction and sealed keypad resist damage from dirt and moisture
- Long 70-hour battery life means fewer recharging cycles
- Capacity tracker prevents overload conditions
- Selectable sample size makes repetitive tasks more efficient
- 2-year warranty

## Applications

- Weighing
- Parts counting
- Percentage weighing
- Dynamic / Animal weighing
- Checkweighing
- Check counting
- Display hold

### About Cole-Parmer:

Cole-Parmer is a leading global manufacturer of proprietary products for fluidics, environmental, analytics, lab essentials, and critical research and process equipment for life sciences. We are proven experts in the fields of fluid handling, environmental testing, cold chain monitoring, sample preparation and handling, and chromatography. Our instruments and components have been used extensively in the biopharma, pharmaceutical, life sciences, environmental, and the food and beverage industries for over 65 years.

Model	IPSW-400-30	IPSW-400-50	IPSW-400-100	IPSW-400-150	IPSW-400-300
Capacity	30lb / 15kg	50lb / 25kg	100lb / 50kg	150lb / 75kg	300lb / 150kg
Readability	0.002lb / 0.001kg	0.002lb / 0.001kg	0.001lb / 0.005kg	0.01lb / 0.005kg	0.02lb / 0.01kg
Pan Size	9.8"x9.8" / 250x250mm		15.7"x11.8" / 400x300mm		23.6"x17.7" / 600x450mm
Weighing Units	g, kg, lb, lb:oz, oz, N, custom unit				
Calibration	External				
Display	Display LCD with 1.6" / 40mm-high digits				
Interface	RS-232				
Power Supply	12VDC 800mA adapter, Internal rechargeable battery				
Operating Temperature	-10° to 40°C / 14° to 104°F				
Housing	Stainless steel				
Overall Dim.	11.6"x13.8"x22.4" / 294x351x568mm		11.8"x20.6"x26.3" / 300x523x667mm		17.7"x27.6"x40.7" / 450x702x1034mm
Net Weight	8kg		11kg		21kg

Supplied by