

# Nordson Divergent™ 3-Way Stopcocks with Bond Port

## Polycarbonate

Cole-Parmer®



Add a stopcock to your system for easy on/off control of your fluid. Three-way stopcocks feature two female Luer fittings plus a bond port. Luer fittings on stopcocks ensure a quick, leak-free connection. The bond port creates a permanent connection; simply match the tube OD to bond port size. Diverter arrow indicates off position, enhancing operational efficiency and minimizing errors.



### CERTIFICATION

ISO 80369-7 compliant

**WARNING:** This product is not approved or intended for, and should not be used for, medical, clinical, surgical, or other patient-oriented applications.

### APPLICATIONS

- Process integration, parallel sampling, and buffer exchange systems
- Ideal for general-purpose, laboratory, and single-use applications

### FEATURES/BENEFITS

- Low-torque handle provides smooth adjustments
- Diverter arrows indicate flow direction
- DEHP free (non-DEHP)
- Transparent polycarbonate body makes it easy to see the flow

USA	+1.800.323.4340
	+1.847.549.7600
Canada	+1.800.363.5900
China	86.21.5109.9909
India	+1.800.266.1244

UK	+44 (0) 1480.272279
	+33 (0) 1.87.170142*
	+49 (0) 937.792030*
Italy	+39.02.84349215
All others	+1.847.549.7600

Cole-Parmer®

[coleparmer.com](http://coleparmer.com)

\*Inquiries from Germany and France are now handled in our St. Neots office by native-speaking experts.

### Specification and Ordering Table

Number of Ports	Port 1	Port 2	Port 3	Bond Port	Max Operating Pressure (psi)	Item Number	Pack Size
3	Female Luer	Female Luer	Bond Port	0.110"	50	<a href="#">56140-93</a>	25
3	Female Luer	Female Luer	Bond Port	0.130"	50	<a href="#">56140-94</a>	25

### FLOW PATH CONFIGURATIONS

