


Filter Adapter Set, Synthetic Rubber



Set of seven synthetic rubber adapter cones designed to fit Buchner and fritted glass funnels, porcelain filters without headpieces, and for reverse flow cleaning of filter crucibles. Adapters can be used to support funnels, Gooch crucibles, filter tubes, and other apparatus in filtration assemblies. These adapters are functional even if not tightly fitted around the funnel stem. Can be used individually or nested with adjacent sizes. Provide an excellent alternative to boring special size holes in rubber stoppers. Autoclavable.

Item No.	Description			
FASN07	Filter Adapter Set of 7, Synthetic Rubber			
Adapter #	Top OD (mm)	Bottom OD (mm)	Top ID (mm)	Bottom ID (mm)
1	21	10	17	6
2	27	16	20	10
3	36	21	30	15
4	45	27	39	21
5	57	37	50	29
6	70	45	62	35
7	90	55	80	46