


Keep your work area clean with wipes from Cole-Parmer



Cole-Parmer® IPA Surface Wipes

General-purpose, disposable 70% isopropyl alcohol (IPA) wipes are ideal for a variety of cleaning applications, including electronics, tools, workspaces, and other surfaces that require quick removal of contaminants. Durable, tear-resistant, lint-free polyester and viscose material is suitable for use in ISO Class 6–7 cleanrooms. Not sterile.

Wipes are presaturated in a solution of 70% isopropyl alcohol (CAS: 63-67-0) and 30% deionized water. The screw-top storage canister includes a convenient flip lid with pull-through dispenser. Each canister contains 100 perforated wipes. Canister is foil-sealed to keep wipes moist before first use.

Wipe size	Quantity	Catalog number
5.5" x 7" (14 x 18 cm)	One canister (100 wipes)	ML-33675-09

Diversey™ Oxivir® Tb Disinfectant Wipes

Tough on pathogens, not on people

- Disinfects surfaces in less than one minute
- Non-irritating to skin and eyes
- Gentle formula suitable for daily use
- Penetrates through the toughest soils
- Kills SARS-CoV-2 (COVID-19 virus)

Oxivir Tb is a one-step disinfectant wipe that is ready to use out of the bottle. Simply wipe any hard, nonporous surface and allow the surface to remain wet for up to one minute. With no added phosphates, phenols, and fragrance, Oxivir Tb cleaner meets the lowest EPA toxicity category and is environmentally responsible—active ingredients break down into oxygen and water.



Disinfectant is effective against key pathogens such as TB, MSRA, Norovirus, and VRE. It kills HBV, HCV, and HIV in one minute; and is a potent tuberculocide and fungicide to effectively clean, disinfect, and deodorize.

Note: Food contact surfaces must be rinsed with potable water after using. Do not use on glassware, utensils, or dishes. It is recommended to test on a small, discreet surface.

Wipe size	Quantity	Catalog number
6" x 7" (15.5 x 17.3 cm)	12 canisters; 160 wipes/ea	ML-78901-04
7" x 8" (17.8 x 20.3 cm)	12 canisters; 60 wipes/ea	ML-78901-06
11" x 12" (28 x 30 cm)	4 canisters; 160 wipes/ea	ML-78901-02

Cole-Parmer® Clean-Wipes™

Won't scratch delicate surfaces

- No residual lint—made of strong, nonwoven cellulose and polyester
- No more struggling with solutions and separate towels—ready to use
- Individually sealed containers for use in cleanrooms

Choose wipes that are premoistened or dry. Convert dry wipes to moistened wipes by adding 8 oz of fluid to container.



Description	Wipe size	Quantity	Catalog number
Premoistened; 70% IPA, 30% DI water	6" x 9" (15.2 x 22.9 cm)	1 canister; 100 wipes/ea	ML-33675-03
Premoistened; 100% DI water			ML-33675-05
Dry			ML-33675-07

Diversey™ Avert® Sporidical Disinfectant Wipes

Cleans, disinfects, and deodorizes all in one step

- Disinfects surfaces in one minute
- Kills *Clostridium difficile* spores in four minutes
- Ready-to-use one-step disinfectant saves you time
- Effective against bacteria, antibiotic resistant bacteria, viruses, and fungi



Avert® is a ready to use, one-step, all-purpose disinfectant cleaner and deodorizer ideal for eliminating *Clostridium difficile* (*C. diff*) spores within four minutes on hard, nonporous surfaces such as glass, glazed porcelain, and stainless steel found in medical and healthcare facilities. The bleach-based formula in these premoistened wipes is also effective within a minute against various pathogens including bacteria, antibiotic resistant bacteria, viruses, and fungi.

Meets minimum active ingredient requirements for CDC, SHEA, APIC, and AORN.

Wipe size	Quantity	Catalog number
6" x 7" (15.5 x 17.3 cm)	12 canisters; 160 wipes/ea	ML-78901-09

Diversey™ Oxivir® Tb Disinfectant Solution

Tough on pathogens, not on people

- Disinfects surfaces in less than one minute
- Ready-to-use one-step disinfectant saves you time
- Non-irritating to skin and eyes

Wipe any hard, nonporous surface and allow the surface to remain wet for up to one minute. With no added phosphates, phenols, and fragrance, this cleaner meets the lowest EPA toxicity category and is environmentally responsible—active ingredients break down into oxygen and water.



Disinfectant is effective against key pathogens such as TB, MSRA, Norovirus, and VRE. It kills HBV, HCV, and HIV in one minute; and is a potent tuberculocide and fungicide to effectively clean, disinfect, and deodorize.

Note: Food contact surfaces must be rinsed with potable water after using. Do not use on glassware, utensils, or dishes. If you are unsure of compatibility, it is recommended to test on a small, discreet surface.

Description	Packaging	Catalog number
Disinfectant solution	Four gallon-sized jugs	ML-78900-87
	One 55-gallon drum	ML-78900-86

ColeParmer.com

USA:
1-800-323-4340
1-847-549-7600

Canada: 1-800-363-5900
China: +86-21-5109-9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: 1-800-266-1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480-272279
All other countries: 1-847-549-7600