

WA-200 Series

Flask Shakers



Eliminate wrist strain

Valuable time can be taken up mixing bottles and flasks by hand. Let our flask shaker take the strain. Holding up to eight flasks or bottles, up to 500ml capacity, it creates a vigorous mixing action by simulating hand shaking. It is especially useful for applications where prolonged shaking is required as it can be left on continuously and won't get tired!

Robust construction, mounted on four rubber feet to absorb vibration and prevent unnecessary movement on the bench. Analogue timer covers 10 to 60 minutes with a manual override. Electronic feed-back control ensures a constant speed irrespective of load.

Supplied with two side arms, eight clamps and Allen key. An optional extension kit is available which consists of two extension arms and four extra clamps increasing the capacity to twelve bottles or flasks.

FEATURES

- Vigorous shaking action simulates hand shaking
- Timed operation or continuous running
- Robust construction
- Ideal for extraction

Specifications

Specifications	WA-200
Speed Range	80 to 800 oscillations/min
Amplitude	1.5 mm
Dimensions (W x D x H)	78 x 27 x 24 cm
Maximum Load	3kg
Net Weight	8.3 kg
Electrical Supply	230 V, 50 Hz, 50 W
IP Rating	31

Ordering Information

Description	Series Number	Model Number	Leagacy SKU	Item Number
Flash shaker with two side arms and eight clamps	WA-200	WA-200	SF1	51900-02
Spare clamp	WA-200	—	SF1/1	99955-87
2 accessory side arms and 4 clamps increases the capacity to twelve bottles or flasks	WA-200	—	SF1/2	99955-88

USA +1.800.323.4340
+1.847.549.7600
Canada +1.800.363.5900
India +1.800.266.1244
China +86.21.5109.9909

UK +44 (0) 1480.272279
+33 (0) 1 87 17 01 42*
+49 (0) 937 792030*
Italy +39.02.84349215
All others +1.847.549.7600

Cole-Parmer®
coleparmer.com

*Enquiries from Germany and France are now handled in our St Neots office by native speaking experts.