Cole-Parmer[®]

WD-200 Series Paraffin Wax Dispenser

- Easy and mess-free tissue embedding
- Heated wax "on tap" for histology and pathology labs
- Enhanced safety with overtemperature cut-out at 120 °C (248 °F)
- Non-drip heated tap to prevent blockages
- Constructed of stainless steel with antimicrobial coating to inhibit bacteria growth
- Slimline design maximizes your bench space
- Outer surface stays cool to the touch
- 3-year warranty



coleparmer.com



The WD-200 series paraffin wax dispenser is designed for on-demand delivery of molten wax. It heats wax to a liquid state and provides convenient 'wax on tap' for histology and pathology laboratories. Large 7.5 L capacity can melt up to 6 kg (13.2 lb) of pelletized wax at a time.

Delivery of molten wax is via a non-drip lever tap with three settings (off, on-demand pull dispense and continuous flow). The inner tank has a filter screen (0.5 mm mesh) fitted to prevent coarse particles from blocking the delivery tap. The wax dispenser tap is heated via a dedicated heating system which prevents blockages and ensures an even flow of wax on-demand.

The wax dispenser has digital temperature control which provides accurate temperature control of the molten wax, and is coupled with an ultra-fast heating system for rapid melting of pelletized wax. The set temperature can be controlled using the up and down arrows on the control panel. The heater indicator will illuminate to show heater activity. When set at 65 $^{\circ}$ C (158 $^{\circ}$ F), the unit takes approximately 90 minutes to fully melt 5 kg (11 lb) of pelletized wax, with significant quantities of molten wax available for use after as little as 30 minutes.

The equipment is protected by a miniature circuit breaker and a temperature cut-out. If an overtemperature condition occurs, the wax dispenser will automatically turn off at 120 $^{\circ}$ C (248 $^{\circ}$ F).



Specifications

Model no.	WD-200		
Capacity	7.5 L		
Material	Stainless steel		
Temperature	Ambient to 70 °C (158 °F)		
Dimensions (W x H x D)	18.1 x 45.5 x 50 cm (7.1 x 17.9 x 19.7 in)		
Warranty	3 years		

Ordering Information

Description	ltem no.	Series no.	Model no.	Legacy SKU
Wax dispenser, 230 VAC, 50/60 Hz	85230-07	WD-200	WD-200D	MH8524
Wax dispenser, 120 VAC, 50/60 Hz	85230-08	WD-200	WD-200D-120	MH8524X1

coleparmer.com

USA: +1.800.323.4340 +1.847.549.7600 Canada: +1 800 363 5900 China: +86 21 5109 9909 France: +33 (0)1 48 63 78 00 Germany: +49 (0) 9377 92030 India: +1 800 266 1244 Italy: +39 (0)2 84349215 UK: +44 (0) 1480 272279 All other countries: +1 847 549 7600

Cole-Parmer[®]