

## Cole-Parmer® Aseptic Sampling Cap For GL45 Glass Media Bottles



The one-piece molded closure system is made with polypropylene cap and TPE tubing. Designed to work with GL45 glass media bottles, the cap provides a closed system to reduce contamination and help protect sample integrity.



### CERTIFICATION

USP Class VI and RoHS compliant

**WARNING:** This product is not approved or intended for, and should not be used for, medical, clinical, surgical, or other patient-oriented applications.

### APPLICATIONS

- Single-use aseptic sampling, tissue and cell culture media, and cell fermentation

### FEATURES/BENEFITS

- Closed system protects sample integrity
- Eliminate cleaning validation
- Reduce assembly and installation time
- Two ports with 1/8" ID x 1/4" OD TPE tubing vent and transfer lines

### Specification and Ordering Table

Cap Material	Cap/Stopper Size	Vent Port Tube Length: in. (cm)	Transfer Port Tube Length: in. (cm)	Bottle Tube Length: in. (cm)	Item Number
Polypropylene	GL45	4" (10)	18" (45)	18" (45)	<a href="#">50142-50</a>

USA	+1.800.323.4340
	+1.847.549.7600
Canada	+1.800.363.5900
China	+86.21.5109.9909
India	+1.800.266.1244

UK	+44 (0) 1480.272279
	+33 (0) 1.87.170142*
	+49 (0) 937.792030*
Italy	+39.02.84349215
All others	+1.847.549.7600

**Cole-Parmer®**  
[coleparmer.com](http://coleparmer.com)

\*Inquiries from Germany and France are now handled in our St. Neots office by native-speaking experts.