

Cole-Parmer® Filling Bell Polycarbonate



Connect tubing to the port at the top of the filling bell and place the filling bell over the receiving vessel for easy aseptic fluid transfer.

The single-use filling bells protect your sterile media against airborne contaminants. The clear polycarbonate filling bells allow for visibility while filling and are lightweight for easy handling.



CERTIFICATION

ASME-BPE, USP Class VI, and RoHS compliant

WARNING: This product is not approved or intended for, and should not be used for, medical, clinical, surgical, or other patient-oriented applications.

APPLICATIONS

- Single-use pharmaceutical, biotech, food, beverage, and other sanitary industries

FEATURES/BENEFITS

- Clear polycarbonate for easy visibility while filling
- Measures 11 cm H x 7.8 cm ID (4-11/32" x 3-1/16") with tubing adapter
- 10 mm OD tube adapter
- Animal-derived ingredient free (ADIF) and animal-derived component free (ADCF)
- Sterilize by autoclave, gamma irradiation, X-ray, or EtO

Specification and Ordering Table

Height: in. (mm)	Material	Color	Outside Diameter (OD): in. (mm)	Inside Diameter (ID): in. (mm)	Item Number
4.3" (110)	Polycarbonate	Clear	3.2" (82)	3.1" (77.6)	50142-95

USA +1.800.323.4340
+1.847.549.7600
Canada +1.800.363.5900
China +86.21.5109.9909
India +1.800.266.1244

UK +44 (0) 1480.272279
+33 (0) 1.87.170142*
+49 (0) 937.792030*
Italy +39.02.84349215
All others +1.847.549.7600

Cole-Parmer®
coleparmer.com

*Inquiries from Germany and France are now handled in our St. Neots office by native-speaking experts.