

METTLER TOLEDO Weighing Paper

Precision in Every Measurement



Effortless Sample Transfer

METTLER TOLEDO weighing paper facilitates easy sample transfer from the balance to the sample container with its smooth surface and flexibility.




Protect from Contamination

Functioning as a clean, disposable layer placed between the sample and the balance's weighing pan, METTLER TOLEDO weighing paper effectively prevents contamination.



Eco-Friendly and Biodegradable

Made from natural fibers, METTLER TOLEDO weighing paper is fully biodegradable, breaking down without leaving harmful residues.



Simplify and Streamline Your Weighing Process Today

Make the most of your weighing workflows with METTLER TOLEDO calendared cellulose weighing paper. Not treated with silicone or any other surfactant, this product is nitrogen-free and made from 100% wood pulp, making it fully biodegradable and food safe.

Weighing paper also reduces the chance of damage or corrosion caused by reactive or sticky materials, preserving the balance's precision and extending its lifespan.

If you are interested in this product, please contact your local METTLER TOLEDO representative or learn more at www.mt.com/lab-accessories.

Material No. 31066128

METTLER TOLEDO