

# Cole-Parmer® Sanitary Clamps

## 304 Stainless Steel

Cole-Parmer®



### Before Choosing a Fixture, Consider the Following

Sanitary fittings are crafted to prevent contamination, ensuring the highest standards of purity and reliability in fluid handling systems. A sanitary fitting consists of three key components: two fittings, a gasket, and a clamp. The gasket is placed between the two fittings, and the clamp holds the assembly together. This design provides a highly sanitary connection that is easy to assemble and disassemble, making it perfect for processes requiring regular maintenance and cleaning.

These easy-to-use clamps are made from 304 stainless steel to provide good chemical resistance.



#### CERTIFICATION

3-A and FDA

### APPLICATIONS

- Great for food, beverage, and dairy applications
- Ideal for cell culture media, aseptic fill and transfer, and sterile filling

### FEATURES/BENEFITS

- Constructed of 304 stainless steel for good chemical resistance
- Compatible with Tri-Clamp® fittings
- Pair with sanitary gasket and fittings for a hygienic connection
- Free technical application support available to assist with fitting and gasket selection

### Specification and Ordering Table

Clamp Size	Item Number	Pack Size
1/2" to 3/4"	<a href="#">50118-42</a>	1
1" to 1-1/2"	<a href="#">50118-43</a>	1
2"	<a href="#">50118-44</a>	1
2-1/2"	<a href="#">50118-45</a>	1
3"	<a href="#">50118-46</a>	1
4"	<a href="#">50118-47</a>	1

USA +1.800.323.4340  
+1.847.549.7600  
Canada +1.800.363.5900  
China +86.21.5109.9909  
India +1.800.266.1244

UK +44 (0) 1480.272279  
+33 (0) 1.87.170142\*  
+49 (0) 937.792030\*  
Italy +39.02.84349215  
All others +1.847.549.7600

Cole-Parmer®

[coleparmer.com](http://coleparmer.com)

\*Inquiries from Germany and France are now handled in our St. Neots office by native-speaking experts.