

Cole-Parmer® Sanitary Gaskets PTFE

Sanitary fittings are crafted to prevent contamination, ensuring the highest standards of purity and reliability in fluid handling systems. A sanitary fitting consists of three key components: two fittings, a gasket, and a clamp. The gasket is placed between the two fittings, and the clamp holds the assembly together. This design provides a highly sanitary connection that is easy to assemble and disassemble, making it perfect for processes requiring regular maintenance and cleaning.

These PTFE gaskets offer excellent chemical, heat, stress, and corrosion resistance for worry-free and leak-free fluid transfer. Gaskets are free of phthalates, TSE/BSE, and animal-derived ingredients, ensuring the integrity of your fluid.



CERTIFICATION

FDA 21 CFR 177.2600, USP Class VI, and 3-A compliant

For additional options and complete offerings, visit [coleparmer.com/c/sanitary-fittings](https://www.coleparmer.com/c/sanitary-fittings)

APPLICATIONS

- Extreme temperature processes and corrosive chemical transfer
- Ideal for general-purpose, laboratory, and single-use applications

FEATURES/BENEFITS

- Excellent chemical and heat resistance for harsh environments
- Operating temperature range of -73 to 260 °C (-100 to 500 °F)
- Compatible with Tri-Clamp® fittings
- Sterilize by autoclave, clean-in-place (CIP), or steam-in-place (SIP)

USA	+1.800.323.4340 +1.847.549.7600
Canada	+1.800.363.5900
China	+86.21.5109.9909
India	+1.800.266.1244

UK	+44 (0) 1480.272279 +33 (0) 1.87.170142*
Italy	+49 (0) 937.792030*
All others	+39.02.84349215 +1.847.549.7600

Cole-Parmer®
[coleparmer.com](https://www.coleparmer.com)

*Inquiries from Germany and France are now handled in our St. Neots office by native-speaking experts.

Specification and Ordering Table

Sanitary Clamp Size	Item Number	Pack Size
1/2"	50133-60	10
3/4"	50133-61	10
1"	50133-62	10
1-1/2"	50133-63	10
2"	50133-64	10
2-1/2"	50133-65	10

PTFE GASKETS

