

# Cole-Parmer®

## Cole-Parmer® Sanitary End Caps Polypropylene



Use sanitary end caps to seal inactive process lines, preventing contamination and leakage. Sanitary end caps are often used in food and beverage, pharmaceutical, and biotechnology applications where maintaining cleanliness is crucial.



### CERTIFICATION

ASME-BPE, USP Class VI, and RoHS compliant

**WARNING:** This product is not approved or intended for, and should not be used for, medical, clinical, surgical, or other patient-oriented applications.

### APPLICATIONS

- Biotech, food and beverage, and pharmaceutical

### FEATURES/BENEFITS

- Cap off fluid system while maintaining cleanliness
- Secure and leak-free sealing
- Animal-derived ingredient free (ADIF) and animal-derived component free (ADCF)
- Polypropylene construction
- Sterilize via gamma irradiation and autoclave

USA	+1.800.323.4340
	+1.847.549.7600
Canada	+1.800.363.5900
China	+86.21.5109.9909
India	+1.800.266.1244

UK	+44 (0) 1480.272279
	+33 (0) 1.87.170142*
	+49 (0) 937.792030*
Italy	+39.02.84349215
All others	+1.847.549.7600

# Cole-Parmer®

[coleparmer.com](http://coleparmer.com)

\*Inquiries from Germany and France are now handled in our St. Neots office by native-speaking experts.

## Specification and Ordering Table

Fitting Material	Sanitary Clamp Size	Item Number
Polycarbonate	3/4"	<a href="#">50142-51</a>
Polypropylene	1-1/2"	<a href="#">50142-53</a>
Polypropylene	2"	<a href="#">50142-55</a>
Polypropylene	3"	<a href="#">50142-57</a>
Polypropylene	4"	<a href="#">50142-59</a>