


Cole-Parmer® Viton® 3-Stop Pump Tubing

Stopped tubing ensures a perfect fit every time. Stops hold tubing in the pump head for optimal tension, resulting in reliable pumping and repeatable flow rates. Three stops give you two tubing sections instead of just one. When one tubing section wears, a new section is ready to use. Viton® tubing features excellent flexibility combined with high chemical resistance. The anti-aging characteristics, expanded temperature ranges, and abrasion resistance offer excellent durability and long service life.

VITON® 3-STOP PUMP TUBING SIZE



Specification and Ordering Table

Dimensions: mm (in.)			Color-Coded Stops	Tubes/ Pack	Item Number
Tubing ID	Length Between Stops	Tubing Length			
1.14 (0.045")	82 (3.2")	400 (15.7")	Red/Red/Red	12	50137-96

APPLICATIONS

- Perfect for fuel and oil lines, aromatic oils, and use in peristaltic pumps
- Ideal for general-purpose, laboratory, and single-use applications

FEATURES/BENEFITS

- Three tubing stops are color-coded for easy size identification
- Sterilization not recommended
- Withstands a wide temperature range of -20 to 200 °C (-4 to 392 °F)
- Flexible tubing with a 60 Shore A durometer
- Opaque black tubing—perfect for sensitive fluids
- Free technical application support available to assist with tubing selection

USA +1.800.323.4340
+1.847.549.7600
Canada +1.800.363.5900
China 86.21.5109.9909
India +1.800.266.1244

UK +44 (0) 1480.272279
+33 (0) 1.87.170142*
+49 (0) 937.792030*
Italy +39.02.84349215
All others +1.847.549.7600

Cole-Parmer®
coleparmer.com

*Inquiries from Germany and France are now handled in our St. Neots office by native-speaking experts.