


Watson-Marlow[®] Peroxide-Cured Silicone 3-Stop Pump Tubing

Peroxide-cured silicone tubing has features that offer optimal durability, with an elongation at break ranging from 440% to 530%, providing a long shelf life and enhanced mechanical strength.

Ideal for use in peristaltic pumps, stops on the tubing are color-coded according to the tubing's inside diameter for easy size identification and ensure a perfect fit on your pump head. Three stops on the tubing offer two pump segments, doubling the life of the tubing.

SILICONE 3-STOP PUMP TUBING



APPLICATIONS

- Perfect for cell culture, reagent delivery, and chromatography systems
- Ideal for general-purpose, laboratory, and single-use applications

FEATURES/BENEFITS

- Designed for use in Watson-Marlow multichannel pump heads
- Specially formulated peroxide-cured silicone offers extended service life
- Exhibits good heat, cold, ozone, UV, and abrasion resistance
- Flexible tubing with a 60 Shore A durometer
- Transparent tubing—perfect for easy flow monitoring
- Free technical application support available to assist with tubing selection

USA	+1.800.323.4340
	+1.847.549.7600
Canada	+1.800.363.5900
China	86.21.5109.9909
French speaker	+33 (0) 1 87 17 01 42

German speaker	+49 (0) 937 792030
India	+1.800.266.1244
Italy	+39.02.84349215
UK	+44 (0) 1480.272279
All others	+1.847.549.7600

Specification and Ordering Table

Brand	Tubing ID: mm	Color-coded stops	Item number	Tubes per pack
Watson-Marlow	0.63	Orange/White	50130-71	6
	0.76	Black/Black	50130-72	6
	0.88	Orange/Orange	50130-73	6
	1.02	White/White	50130-74	6
	1.14	Red/Red	50130-75	6
	1.29	Gray/Gray	50130-76	6
	1.42	Yellow/Yellow	50130-77	6
	1.52	Yellow/Blue	50130-78	6
	1.65	Blue/Blue	50130-79	6
	1.85	Green/Green	50130-80	6
	2.05	Purple/Purple	50130-81	6
	2.29	Purple/Black	50130-82	6
	2.54	Purple/Orange	50130-83	6
	2.79	Purple/White	50130-84	6



The information contained in this document is believed to be correct but Cole-Parmer accepts no liability for any errors it contains, and reserves the right to alter specifications without notice. Watson-Marlow is a registered trademark of Watson-Marlow Limited. Please state the item number when ordering.

Cole-Parmer[®]
an antylia scientific company