

Cole-Parmer® Tube Rotators

TR-200 and TR-200D



Rotating shakers provide smooth effective motion for your samples..

These tube rotators feature an adjustable angle of rotation from horizontal for minimal mixing to vertical for full end-over-end mixing. An integrated drip tray prevents workspace spills from the rotating tubes and both units offer a choice of six tube holders for most applications. The tube holders can be removed and replaced quickly and easily. These tube rotators can use two tube holders simultaneously with our dual holder accessory (00077-23). The rotators can be used in incubators up to 60 °C and in cold rooms down to 4 °C.

TR-200

- Gives gentle but effective mixing, ideal for keeping biological samples in suspension.
- The rotator has a constant speed of 20 rpm.

TR-200D

- Ideal for aerating cultures, keeping biological samples in suspension and for general mixing applications.
- The rotator has a variable speed of 2 to 40 rpm for gentle rolling or vigorous mixing of samples.
- A digital timer and speed display allows procedures to be accurately repeated for optimal results.
- Scale to measure mixing angle for future reference or continuity.

FEATURES

- Adjustable angle of rotation mixes samples effectively with gentle or vigorous motion
- A variety of tube holder discs allows customization for your application
- Integrated drip tray prevents workspace spills



TR-200D Tube Rotator

Specifications

Specification	TR-200	TR-200-120	TR-200D	TR-200D-120
Width (cm)	22.5	22.5	22.5	22.5
Height (cm)	29.0	29.0	29.0	29.0
Depth (cm)	31.0	31.0	31.0	31.0
Weight (kg)	5.2	5.2	5.3	5.3

Ordering Information

Description	Series Number	Model Number	Legacy Cole-Parmer SKU	Legacy Stuart SKU	Item Number
Rotator, analog control, fixed speed, 230 VAC, 50 Hz	TR-200	TR-200	07650-11	SB2	51901-26
Rotator, analog control, fixed speed, 120 VAC, 60 Hz	TR-200	TR-200-120	07650-07	SB2/120V/60	51901-27
Rotator, digital control, variable speed, 230 VAC, 50 Hz	TR-200	TR-200D	07650-13	SB3	51901-28
Rotator, digital control, variable speed, 120 VAC, 60 Hz	TR-200	TR-200D-120	07650-17	SB3/120V/60	51901-29

USA +1.800.323.4340
 +1.847.549.7600
 Canada +1.800.363.5900
 India +1.800.266.1244
 China 86.21.5109.9909

UK +44 (0) 1480.272279
 +33 (0) 1 87 17 01 42*
 +49 (0) 937 792030*
 Italy +39.02.84349215
 All others +1.847.549.7600

Cole-Parmer®
 coleparmer.com

*Enquiries from Germany and France are now handled in our St Neots office by native speaking experts.