



# Cole-Parmer®

## Cleanroom Delicate Task Wipers

PROTECT CRITICAL GLASS SURFACES, LENSES AND INSTRUMENTS FROM CLEANING-INDUCED SCRATCHES AND LINT



### KIMTECH Science\* Kimwipes\* EX-L

These Clean-room Wipers from Cole-Parmer are the first choice for Scientists and Professionals for maintaining critical glassware and instruments. The no-lint, single-ply wipers are produced from proprietary non-woven, finely textured tissue that is 2.8 times more absorbent than commercially available paper towels – reducing consumption and disposal costs. Will not scratch delicate surfaces even when used wet. Supplied in unique pop-up boxes with antistatic Lintguard® polyshield that further assist in eliminating lint and static buildup. Available in cases for best value.

#### Specifications :

- Manufactured from 100% virgin fiber
- Multi-fold dispensing format enables one-at-a-time dispensing
- Compliance : Class 100
- Available in two sizes

Catalog Number	Description	Wipe Size	Qty/bx	Box/cs	Price/cs
YK-33670-04	Kimwipes EX-L, 1 ply	4.5" X 8.5"	280	60	Rs. 13,500
YK-33670-06		15" X 16.75"	140	15	Rs. 11,000

*Prices are Ex-Warehouse, Taxes Extra, Valid upto Dec 2008*

\* Trademark of Kimberly-Clark Worldwide, Inc.