

CERTIFICATE OF QUALITY

This document certifies that the bellow product meets the following criteria:

| | |
|----------------------|---|
| Catalog Number: | 06027-56 |
| Product Description: | Amber Wide-Mouth Boston Round Bottle-HDPE 250ml |
| Lot Number: | XXXXXXXXXX |
| Date of Manufacture: | XXXXXXXXXX |
| Country Of Origin: | India. |

Quality System Compliance: Products are manufactured, inspected and found to be in compliance with product and quality specification requirements as documented in the ISO 9001:2015 & ISO13485:2016 compliant Quality Management System. Inspection records are reviewed and signed for product release.

QC Testing: Representative production samples are collected and inspected in accordance with current applicable standard protocols.

| Test Name | Test Reference | Accuracy | Result |
|----------------|----------------|----------------------------------|------------------|
| 1) Leak Test | QWI-QC-03 | General Insp. Level- I, AQL- 1.5 | Leak proof |
| 2) Autoclaving | QWI-QC-07 | General Insp. Level- I, AQL- 1.5 | NA |
| 3) Dimensional | QSF-QCIR-02 | General Insp. Level- I, AQL- 1.5 | Within tolerance |

Product Specifications:

Material: Bottle- High Density Polyethylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.
Closure- Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.

BSE/TSE Statement: The above mentioned product is manufactured without the use of Raw materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge.

Non-Pyrogenic: Not applicable.

Dnase and Rnase Free: Not applicable.

Cytotoxicity: Not applicable.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Issued on XXXXXXXXX.