

CERTIFICATE OF QUALITY

This document certifies that the bellow product meets the following criteria:

Catalog Number:	06036-01, 06036-03
Product Description:	Economy Amber Wide-Mouth Boston Round Bottle-HDPE 30ml
Lot Number:	XXXXXXXXXX
Date of Manufacture:	XXXXXXXXXX
Country Of Origin:	India.

Quality System Compliance: Products are manufactured, inspected and found to be in compliance with product and quality specification requirements as documented in the ISO 9001:2015 & ISO13485:2016 compliant Quality Management System. Inspection records are reviewed and signed for product release.

QC Testing: Representative production samples are collected and inspected in accordance with current applicable standard protocols.

Test Name	Test Reference	Accuracy	Result
1) Leak Test	QWI-QC-03	General Insp. Level- I, AQL- 1.5	Leak proof
2) Autoclaving	QWI-QC-07	General Insp. Level- I, AQL- 1.5	NA
3) Dimensional	QSF-QCIR-02	General Insp. Level- I, AQL- 1.5	Within tolerance

Product Specifications:

Material: Bottle- High Density Polyethylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.
Closure- Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.

BSE/TSE Statement: The above mentioned product is manufactured without the use of Raw materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge.

Non-Pyrogenic: Not applicable.

Dnase and Rnase Free: Not applicable.

Cytotoxicity: Not applicable.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Issued on XXXXXXXXX.