

Part Number	06264 - 22
Description	Micro Spatula VN, Sterile, PS, .3mL; 100/PK
Lot Number	104702
Expiry Date	25 April 2027
Material	Polystyrene (PS)
Certificate Issued	14 September 2022

- We confirm that the goods are manufactured with Polystyrene that conforms to:
 - FDA CFR 177.1640
 - EU regulations 10/2011
 - EC Regulations 1935/2004
- The Polystyrene resin used in the manufacture of the above products is BSE/TSE free. During manufacture this item does not come into contact with any products or processes that may pose a BSE/TSE risk.
- These goods are moulded and packed in a class 100,000 cleanroom. They are individually packed.
- This product is designed for single use only. It must not be used multiple times.
- The information detailed above has been given to us by our suppliers, it is passed on in good faith and is to the best of our knowledge correct.