

<b>Part Number</b>	06286 - 21
<b>Description</b>	Slot Sampler, 1 opening, Sterile, HDPE, 620mm length, 50mL; 20/PK
<b>Lot Number</b>	104709
<b>Expiry Date</b>	27 May 2027
<b>Material</b>	HDPE (High Density Polyethylene)
<b>Certificate Issued</b>	14 September 2022

---

- We confirm that the goods are manufactured with HDPE that conforms to:
  - FDA CFR 177.1520
  - EU regulations 10/2011
  - EC regulations 1935/2004
- The HDPE resin used in the manufacture of the above products is BSE/TSE free. During manufacture this item does not come into contact with any products or processes that may pose a BSE/TSE risk.
- These goods are assembled and packed in a class 100,000 cleanroom. They are individually packed.
- This product is designed for single use only. It must not be used multiple times.
- Ultrasonic welding is used to adhere the components
- The information detailed above has been given to us by our suppliers, it is passed on in good faith and is to the best of our knowledge correct.