

Part Number	06286-30
Description	Powder Sample Collector, Sterile, HDPE, 500mm length, 10mL; 20/PK
Lot Number	104711
Expiry Date	20 May 2027
Material	HDPE (High Density Polyethylene)
Certificate Issued	14 September 2022

- We confirm that the goods are manufactured from HDPE that conforms to:
 - FDA CFR 177.1520
 - EU Regulations 10/2011
 - EC Regulations 1935/2004
- The HDPE resin used in the manufacture of the above products is BSE/TSE free. During manufacture this item does not come into contact with any products or processes that may pose a BSE/TSE risk.
- These goods are assembled and packed in a class 100,000 cleanroom. They are individually packed.
- This product is designed for single use only. It must not be used multiple times.
- Ultrasonic welding is used to adhere the components

The information detailed above has been given to us by our suppliers, it is passed on in good faith and is to the best of our knowledge correct.