

## CERTIFICATE OF QUALITY

**Description** : Media Bottle, PET, Sterile, 1000 mL  
**Catalog Number** : 06287-03  
**Material** : Cap: Polyethylene / Flask: Polyethylene terephthalate  
**Lot Number** : BPETS22090101  
**Expiry Date** : 2025-06  
**Visual appearance** : Within acceptability limits  
**Physical Properties** :

**Sterilization:** Products had been sterilized by Gamma Irradiation Sterilization  
Assurance Level: SAL 10<sup>-6</sup> Irradiation dosage complied with ISO 11137.

**DNase/RNase Free:** Products are tested to be free of DNase/RNase contamination.

**Non-pyrogenic:** Products are tested to meet specification of Non-pyrogenic,  
less than 0.5 Eu/ml.

### *Other Information*

**Appearance:** Pass

**Dimensions:** Pass

**Packaging:** Pass