

## TECHNICAL DATA SHEET

<b>Description</b>	: Media Bottle, PET, Sterile, 1000 mL
<b>Catalog Number</b>	: 06287-03
<b>Volume</b>	: 1000 ml
<b>Neck Size</b>	: 32 mm
<b>Dimensions</b>	: 93 x 93 x 214 mm
<b>Material</b>	: Cap: Polyethylene / Flask: Polyethylene terephthalate
<b>Color</b>	: Transparent
<b>Apperance</b>	: Pass
<b>Dimensional Control</b>	: Pass