



 Cole-Parmer®

TERRASPRAY™ RINSE STATIONS

INCREASE YOUR PRODUCTIVITY WHILE MINIMIZING WASTE



- Compact design keeps water contained—no more wet floors and slippery hazards!
- Choose from countertop or floor models
- Avoid expensive automated systems
- Clean all sizes and styles of vessels and carboys up to 50 gallons (200 liters)

These rinse stations eliminate manual cleaning of vessels and carboys, allowing you to clean and sanitize quickly and without fuss. Easy-to-use systems set up quickly and do not require electricity. Rinse stations connect directly to water inlet and drain for easy washing and cleanup.

Countertop unit handles vessels with a diameter up to 15" (38.1 cm) and features a hand-controlled valve on the side of the unit.

Floor model accommodates vessels with a maximum diameter of 19" (48.3 cm) and has a foot-controlled valve to turn water on and off.



**FREE TECHNICAL
APPLICATIONS ASSISTANCE!**

In the US

625 E. Bunker Court
Vernon Hills, IL 60061
Call toll-free 800-323-4340
Phone: 847-549-7600
Fax: 847-247-2929
www.coleparmer.com

In Canada

Call toll-free 800-363-5900
Phone: 514-355-6100
Fax: 514-355-7119
www.coleparmer.ca

In the United Kingdom

Free phone: 0500-345-300
Phone: 020-8574-7556
Fax: 020-8574-7543
www.coleparmer.co.uk

In India

Phone: 91-22-6716-2222
Fax: 91-22-6716-2211
www.coleparmer.in

International customers

Call 847-549-7600 to reach our
International Sales Department
or contact your local dealer.

PERFORMANCE SPECIFICATIONS

Materials

Top, base, and legs: polyethylene
Valves: PP, Viton® and Kynar® PVDF coated stainless steel springs
Spray head: PVDF

Maximum vessel diameter

Countertop model: 15" (38.1 cm)
Floor model: 19" (18.3 cm)

Rinse control

Countertop model: hand-controlled valve
Floor model: foot-controlled valve

Water supply inlet: ½" NPT(F)

Drain line: 1½" NPT(F)

Catalog number	Washer type	Dimensions (W x H x D)
ML-98965-06	Countertop	19" x 27½" x 19" (48.3.3 x 70 x 48 cm)
ML-98965-08	Floor	22" x 53" x 22" (55.9 x 134.6 x 55.9 cm)

**For further information,
go to ColeParmer.com
or call 800-323-4340.**