

Transport critical samples with ease and less worry



- Improve sample protection—maintain specific temperature ranges for extended time periods
- Choose from durable transport bags or boxes to meet your needs
- Modular, color-coded cooling packs allow flexible configurations for samples
- Maintain sample integrity with optional plastic security seals sold separately

PolarSafe® Transport Bags

- Economical, lightweight insulated nylon bag for hand carrying samples
- Shoulder strap for easier handling
- Keep sample documents safe with integrated pocket
- Sizes: 5 L, 10 L, 16 L, and 30 L



PolarSafe Transport Boxes

- Extremely durable, polypropylene insulated walls for improved insulation
- Space-saving stackable design
- Move boxes easily by using integrated handles or shoulder strap
- Sizes available: 5 L, 10 L, and 20 L



PolarSafe Modular Cooling Packs

- Maintain temperature ranges for extended periods of time
- Cooling end-caps, cooling frames, and cooling blocks are stackable and fit all internal bag/box volumes
- Available in +4°C and +22°C temperatures
 - The 4°C packs maintain 2 to 8°C for 8 to 15 hours*
 - The 22°C packs maintain 20 to 24°C for 8 to 15 hours*
- Color coded for convenience—4°C packs are blue, 22°C packs are green



Cooling end-caps



Cooling frames



Cooling block

*Based on 25°C ambient temperature. Estimates can vary depending on how samples are packed. Additional temperature/time data available at ColeParmer.com

PolarSafe® Transport Kits

- Transport bags and boxes are conveniently bundled with recommended cooling packs, saving you money while protecting your samples



Ordering Information

Transport Bags

Description	Catalog number
5-L capacity	ML-03726-00
10-L capacity	ML-03726-01
16-L capacity	ML-03726-02
30-L capacity	ML-03726-03

Transport Boxes

Description	Catalog number
5-L capacity	ML-03726-04
10-L capacity	ML-03726-05
20-L capacity	ML-03726-06

Cooling Packs

Description	Catalog number
End-caps, 22°C, green. Pack of 2	ML-03726-07
Frames, 22°C, green. Pack of 2	ML-03726-08
End-caps, 4°C, blue. Pack of 2	ML-03726-09
Frames, 4°C, blue. Pack of 2	ML-03726-12
Block, 22°C, green, 1 L	ML-03726-13
Block, 22°C, green, 1 L, slim	ML-03726-14
Block, 22°C, green, 2 L	ML-03726-18
Block, 4°C, blue, 1 L	ML-03726-15
Block, 4°C, blue, 1 L, slim	ML-03726-16
Block, 4°C, blue, 2 L	ML-03726-17

Security Seals

Description	Catalog number
Plastic security sealing ties with blank label area. Case of 100	ML-03726-19

Transport Bag Kits

Description	Catalog number
5-L bag with two 4°C end-caps and two 4°C frames	ML-03726-32
5-L bag with two 22°C end-caps and two 22°C frames	ML-03726-33
10-L bag with two 4°C end-caps and four 4°C frames	ML-03726-34
10-L bag with two 22°C end-caps and four 22°C frames	ML-03726-35
16-L bag with four 4°C blocks (1 L)	ML-03726-36
16-L bag with four 22°C blocks (1 L)	ML-03726-37
30-L bag with two 4°C blocks (1 L, slim) and four 4°C blocks (1 L)	ML-03726-38
30-L bag with two 22°C blocks (1 L, slim) and four 22°C blocks (1 L)	ML-03726-39

Transport Box Kits

Description	Catalog number
5-L box with two 4°C end-caps and two 4°C frames	ML-03726-40
5-L box with two 22°C end-caps and two 22°C frames	ML-03726-41
10-L box with two 4°C end-caps and four 4°C frames	ML-03726-42
10-L box with two 22°C end-caps and four 22°C frames	ML-03726-43
20-L box with two 4°C blocks (2 L)	ML-03726-44
20-L box with two 22°C blocks (2 L)	ML-03726-45

Note: These bags and boxes can be used as outer packaging according to P650 per UN3373 (category B).*

***WARNING** – product is suitable as outer packaging only according to P650 per UN3373 (category B).

To comply with shipping regulations for dangerous goods and infectious substances, additional packaging materials will be needed and packages containing infectious substances must be labelled in compliance with the U.S. Department of Transportation's Hazardous Materials Regulations and the International Air Transport Association Dangerous Goods Regulations and other applicable laws and regulations. Failure to comply with these requirements may be violation of federal law.

The use of this product to transport dangerous goods or infectious substances should only be done by a certified dangerous goods specialist.

ColeParmer.com

USA:
1-800-323-4340
1-847-549-7600

Canada: 1-800-363-5900
China: +86-21-5109-9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: 1-800-266-1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480-272279
All other countries: 1-847-549-7600