

## CERTIFICATE OF QUALITY

This document certifies that the bellow product meets the following criteria:

Catalog Number:	62300-47
Product Description:	WM Wash Bottle LDPE Safety Labeled Vented -Distilled Water 1000ml
Lot Number:	XXXXXXXX
Date of Manufacture:	XXXXXXXX
Country Of Origin:	India.

Quality System Compliance: Products are manufactured, inspected and found to be in compliance with product and quality specification requirements as documented in the ISO 9001:2015 & ISO13485:2016 compliant Quality Management System. Inspection records are reviewed and signed for product release.

**QC Testing:** Representative production samples are collected and inspected in accordance with current applicable standard protocols.

Test Name	Test Reference	Accuracy	Result
1) Leak Test	QWI-QC-03	General Insp. Level- I, AQL- 1.5	Leak proof
2) Autoclaving	QWI-QC-07	General Insp. Level- I, AQL- 1.5	NA
3) Dimensional	QSF-QCIR-03	General Insp. Level- I, AQL- 1.5	Within tolerance

### Product Specifications:

**Material:** Bottle- Low Density Polyethylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.  
Closure- Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.

**BSE/TSE Statement:** The above mentioned product is manufactured without the use of Raw materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge.

**Non-Pyrogenic:** Not applicable.

**Dnase and Rnase Free:** Not applicable.

**Cytotoxicity:** Not applicable.

**Latex Free Statement:** No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Issued on XXXXXXXX.