

Support Plates

Round Support Plate

Aluminum plate supports beakers, mantles and petri dishes. Fits standard lab-frames using Talboys connectors or holders, shown on pages 14-16.

Plate Diameter	Distance from Rod	Part Number
6" (152mm)	2.6" (66mm)	08027-16



Support Plates

Designed to hold Hotplates/Stirrers and other apparatus to lab-frames or ring stands. Aluminum construction offers strength and durability. Support plates include a non-skid rubber mat and a built-in holder that grips rods up to 0.75" (19mm) in diameter.

Size	Plate Dimensions	Distance from Rod	Part Number
Small	9.5 x 6.5" (241 x 165mm)	1.9" (48mm)	08059-10
Medium	12.75 x 9" (324 x 229mm)	1.9" (48mm)	08059-20
Large	15.75 x 12" (400 x 305mm)	1.9" (48mm)	08059-30



In The Know: Get Free Technical Assistance

Choose the right product for your application

Contact an Application Specialist

- Call 800-323-4340
- e-mail techinfo@coleparmer.com
- Visit ColeParmer.com/tapps
- Live Ch@t Online Support: Connect directly for immediate answers