

CLAMPS & SUPPORTS

Ultra Flex Support Systems

Ultra Flex Support Systems have unique flexible arms that are extremely versatile and can be placed in virtually any position or angle. The Ultra Flex comes in three different systems; base plate, lab-frame connector, and bench clamp. Each system includes a 2-Prong Clamp Head, 3-Prong Clamp Head, Spring Clamp Head, nickel-plated flex arm in 305mm (12") or 457mm (18") length, and utility wrench. 2-prong and 3-prong clamp heads also supplied with non-slip vinyl sleeves, and for temperatures above 100°C, fiberglass covers.

Arm Length,	Ultra Flex 12:	305mm (12")
	Ultra Flex 18:	457mm (18")
Arm Diameter:		13mm (0.51")
Max. Grip Size,	2-Prong Clamp Head:	75mm (2.95")
	3-Prong Clamp Head:	69mm (2.72")
	Spring Clamp Head:	13mm (0.51")

2-Prong Clamp Head



3-Prong Clamp Head



Spring Clamp Head



Ultra Flex Support System with Base Plate

- **Stable base plate**
- **All metal construction**

Ultra Flex Support System with Base Plate is designed with a stable, painted steel base plate with chemical resistant black finish, that fits easily on bench tops or in fume hoods. Base plate measures 127 x 127 x 12.7mm (5 x 5 x 0.5").

Description	Part Number
Ultra Flex 12 with Base Plate	08029-00
Ultra Flex 18 with Base Plate	08029-02



08029-02

Ultra Flex Support System with Lab-Frame Connector

- **Cast alloy lab-frame connector**
- **Mounts to lab-frames or support stands**

Ultra Flex Support System with Lab-Frame Connector securely attaches to support stands, lab-frames, or any support rod up to 19mm (0.75") in diameter. Ideal for use in fume hoods. Lab-frame connector is made of cast alloy.

Description	Part Number
Ultra Flex 12 with Lab-Frame Connector	08029-04
Ultra Flex 18 with Lab-Frame Connector	08029-06



08029-04

Ultra Flex Support System with Bench Clamp

- **Aluminum bench clamp**
- **Easily attaches to most lab benches**

Ultra Flex Support System with Bench Clamp easily attaches to the side of your lab bench or

Description	Part Number
Ultra Flex 12 with Bench Clamp	08126-36
Ultra Flex 18 with Bench Clamp	08126-37



08126-37

Ultra Flex Support System Replacement Parts

Description	Part Number
2- Prong Clamp Head	08126-39
3- Prong Clamp Head	08126-44
Spring Clamp Head	08126-46
12" Flex Arm	08126-38
18" Flex Arm	08129-43
Base Plate	08129-40
Lab-Frame Connector	08129-42
Bench Clamp	08129-41