

Quality and Accuracy

Thermo Scientific Orion ROSS pH Electrodes



The best choice for your pH measurement needs, Thermo Scientific™ Orion™ ROSS™ pH electrodes offer unsurpassed accuracy and

trouble-free measurements in a variety of sizes and styles for every application.

All ROSS pH electrodes:

- deliver improved measurements
- are compatible with any brand of meter
- are TRIS compatible for use with most samples, including proteins and sulfides
- are mercury and silver free to prevent sample contamination

ROSS Ultra™ pH electrodes offer all the advantages of ROSS electrode technology plus a two year replacement warranty – double the warranty period of most electrodes.

ROSS pH electrodes offer the best accuracy, stability and response found in a pH electrode.

ROSS Ultra pH/ATC Triode™ electrodes include a built-in temperature sensor to measure both pH and temperature with one electrode, minimizing the amount of required sample.

ROSS Sure-Flow™ pH electrodes have a unique free flowing junction that prevents the electrode from clogging and greatly simplifies electrode cleaning and maintenance.

ROSS Micro pH electrodes utilize a 3mm tip diameter to measure samples as small as 15 µL and is more durable than standard micro pH electrodes.

Free!

ROSS Storage Solution is specially formulated to optimize performance of the ROSS pH electrode and prolong its life. (Part No. 810001)

Orion Storage Sleeve and Base protects your ROSS pH electrode from breakage when not in use and it won't tip over due to its weighted removable base. As an additional convenience, the pH bulb stays conditioned while in storage so that it's ready to use when removed from the sleeve. (Part No. 810017)

Purchase a new ROSS pH electrode and receive 2 bottles of ROSS storage solution and an electrode storage sleeve for free!

Thermo
SCIENTIFIC

Part No.	Description
ROSS Ultra pH Electrodes	
8102BNUWP	ROSS Ultra Glass Refillable pH Electrode
8103BNUWP	ROSS Ultra Semi-Micro Glass Refillable pH Electrode
8104BNUWP	ROSS Ultra Rugged Bulb Glass Refillable pH Electrode
8115BNUWP	ROSS Ultra Semi-Micro Epoxy Refillable pH Electrode
8135BNUWP	ROSS Ultra Flat Surface Epoxy Refillable pH Electrode
8156BNUWP	ROSS Ultra Epoxy Refillable pH Electrode
ROSS Ultra pH-ATC Triode Electrodes	
8107BNUMD	ROSS Ultra Triode Epoxy Low Maintenance Gel pH-ATC Electrode
8107UWMMD	ROSS Ultra Triode Epoxy Low Maintenance Gel pH-ATC Electrode, 3 Meter Cable
8157BNUMD	ROSS Ultra Triode Epoxy Refillable pH-ATC Electrode
8157UWMMD	ROSS Ultra Triode Epoxy Refillable pH-ATC Electrode, 3 Meter Cable
8302BNUMD	ROSS Ultra Triode Glass Refillable pH-ATC Electrode
ROSS Sure-Flow pH Electrodes	
8165BNWP	ROSS Sure-Flow Epoxy Refillable pH Electrode
8172BNWP	ROSS Sure-Flow Glass Refillable pH Electrode
8175BNWP	ROSS Sure-Flow Semi-Micro Epoxy Refillable pH Electrode
ROSS Micro pH Electrode	
8220BNWP	ROSS Micro Glass Refillable pH Electrode
ROSS pH Electrodes	
8102BN	ROSS Glass Refillable pH Electrode
8103BN	ROSS Semi-Micro Glass Refillable pH Electrode
8104BN	ROSS Rugged Bulb Glass Refillable pH Electrode
8115BN	ROSS Semi-Micro Epoxy Refillable pH Electrode
8135BN	ROSS Flat Surface Epoxy Refillable pH Electrode
815600	ROSS Epoxy Refillable pH Electrode
8163BNWP	ROSS Spear Tip Glass Refillable pH Electrode



Take advantage of this limited-time ROSS pH Electrode Promotion

Offer valid April 1, 2013 through September 30, 2013

1. Purchase a qualifying Thermo Scientific Orion ROSS pH electrode.
2. Complete the form below and save the file to your computer.
3. Return the form: Email to **wateranalysis@thermofisher.com** along with a scanned copy of the original invoice or packing slip or Fax to (978) 215-4013 along with the invoice/packing slip.

Contact Name: _____

Company Name: _____

Shipping Address: _____

Telephone No: _____

E-mail Address: _____

I would like to receive e-mail updates about the latest Thermo Scientific Orion products: Yes No

What brand are you currently using?

For your pH meters: _____


For your pH electrodes: _____

For your pH buffers: _____

What's your application? _____

Remember to include your invoice or packing slip

For more details regarding this promotion contact your local sales representative or call us at 1-978-232-6000 or 1-800-225-1480. The offer applies to new qualifying ROSS pH electrodes purchased between April 1, 2013 and September 30, 2013. Redemption form must be submitted by October 31, 2013. This offer cannot be combined with any other offer or discount. Limit five promotion redemptions per customer, company and/or shipping address. Please allow 6-8 weeks for delivery of promotional items. This offer is void where prohibited by law or company policy.



thermoscientific.com/water

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Water Analysis Instruments

North America
Toll Free: 1-800-225-1480
Tel: 1-978-232-6000

Netherlands
Tel: (31) 033-2463887
info.water.uk@thermo.com

India
Tel: (91) 22-4157-8800
wai.asia@thermofisher.com

Japan
Tel: (81) 045-453-9175
wai.asia@thermofisher.com

China
Tel: (86) 21-68654588
wai.asia@thermofisher.com

Singapore
Tel: (65) 6778-6876
wai.asia@thermofisher.com

Australia
Tel: (613) 9757-4300
in Australia (1300) 735-295
InfoWaterAU@thermofisher.com

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific